

novaRAY



Advanced interventional imaging

Cardiac fluoroscopy




X-ray angiogram
Diagnostic gold standard


Cardiac fluoroscopy



X-ray angiogram
Diagnostic gold standard



Balloon angioplasty




Stents

Guide all interventions

novaRAY

\$1.8B interventional imaging market

Cardiac
Catheterization Systems
\$1.1B



Neurovascular
Catheterization Systems
\$0.1B

Abdominal & Peripheral Vascular
Catheterization Systems
\$0.6B

Cardiac Catheterization Systems


□ Average price: \$1M

□ Service life: 8 years

Source: Frost & Sullivan

novaRAY

Conventional fluoroscopy system



- Point x-ray source
- Large-area detector
- Shadowgram image

Conventional fluoroscopy limitations



- ❑ Excessive radiation
- ❑ Poor image quality in large patients / steep angles
- ❑ Anatomical clutter
- ❑ Cumbersome quantitative vessel sizing
- ❑ Obstructed access to patient

NovaRay solution

Entirely new
fluoroscopy technology

- Dramatically reduced radiation
- Enhanced imaging
 - Continuous diagnostic quality
 - Superior quality in large patients
- Unique real-time 3D multi-slice tomography
 - Enhanced visualization
 - Simple accurate vessel sizing
- Open patient access



novaRAY

NovaRay fluoroscopy technology


New architecture
and subsystems

- ▮ Inverse geometry
- ▮ Scanning x-ray source
- ▮ High-efficiency detector array
- ▮ Image-reconstruction computer



novaRAY


NovaRay system




- ❑ Inverse geometry
- ❑ Scanning x-ray source
- ❑ Small-area detector
- ❑ Collect many small overlapping images
- ❑ Scan entire field of view 30 times / second
- ❑ Reconstruct full images in real time

NovaRay proprietary imaging chain


High-efficiency photon-counting x-ray
detector arrays




High-power
scanning x-ray source



High-speed multi-slice tomosynthesis
reconstruction




NovaRay technology





novaRAY

Scanning x-ray source



- 25 kW continuous operation
- 70 kVp - 120 kVp
- 0.4 mm beam spot
- $\sim 100 \text{ kW/mm}^2$
- 23 cm x 23 cm
- 100 x 100 collimator holes
- 1 μs dwell, 0.25 μs move time
- 30 frames per second
- No moving parts in vacuum
- No output decline over time


Photon-counting detector



- Solid state, room temperature detector
- Photon counting
 5×10^6 photons / sec / mm²
- 1.25 μ s cycle time
- High-efficiency 90% DQE(0)

Photon-counting detector array

4 x 8 array of detector hybrids



←→ 1.35 cm


- 5.4 x 10.8 cm active area
- 0.33 mm pixel
- 160 x 320 (51,200) detector elements
- $100 \times 100 \times 160 \times 320 = 5 \times 10^8$ "rays"

Image-Reconstruction Requirements

- Back projection reconstruction algorithm
- 40 Gbit / sec input from 640 x 320 sensor array
 - Condensed to 10 Gbps for our case
- 1 Gpixel / sec image reconstruction
 - Multiple 1,000 x 1,000-pixel focal planes at 30 fps
- 1M pixel “best-focus” image (SIFT) at 30 fps
- ~1.4 TOPS (for 32 planes)

Low x-ray scatter

Conventional



NovaRay



Higher signal-to-noise ratio = improved image quality

novaRAY

Exposure reduction factors

- ▣ Five factors contribute to exposure reduction
 - Low x-ray scatter
 - No anti-scatter grid
 - Higher detector efficiency
 - Adaptive scanning (new)
 - Larger entrance area

Unique multi-slice tomography

Conventional
Shadowgram



NovaRay
Reconstruct multiple slices




novaRAY

Single-slice imaging




NovaRay



novaRAY

Volume imaging



Anatomical clutter removed

NovaRay



Selectable slices

novaRAY